FABRICANTES FELIGNEO, S.A. DE C.V.  
Alfonso Garzón Santibañez No. 7 , Col. Indígena San Juan de Ocotán  
C.P. 45019, Zapopan, Jalisco, México.  
Tel. 01 33) 31 10 60 02, 31 10 60 77, Fax. 31 10 61 03  
[Http://www.felisa.com.mx](http://www.felisa.com.mx)

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## OPERATION INSTRUCTIONS FOR STERILIZERS

**IMPORTANT:** The care you take in reading and following this instructions will determinate the satisfactory service you will receive from your Sterilizer. This equipment must be utilized only for sterilization processes, any other application will invalidate the guarantee.

**UNPACKING AND INSPECTION.** Carefully remove the sterilizer from the shipping case and check it, this equipment is specially packed with polyurethane foam, but if damage has occurred from shipment a claim must be filled with the carrier immediately; preserve the shipping container for inspection by the carrier and contact your dealer or Felisa.

**INSTALLATION.** Read the specifications tag in the cabinet for check the required voltage and potency and plug-in your Sterilizer to the adequate socket. We recommend that this socket be independent from others and be supplied from the principal electric line. The Sterilizer must be installed on the floor, leaving a minimum space of 20 cm between the Sterilizer and any other vertical surface, keep the area around free of any material to allow ventilation in the area of the bottom. **CAUTION:** In case that steam will be released inside the place, steer the steam exit to a zone where it cannot damage somebody or something.

**INTRODUCTION.** All FELISA Sterilizers have two security devices, former a security valve (18) on the evaporator lid (this valve can liberate excessive pressure) and latter a pressure switch (38) that cut off the electric power in heating elements, both are calibrated in our plant and will act always that the safe pressure limit is exceeded. The ANALOGIC MODELS have a knob (11) for the temperature control and the user himself must operate the steam valve and count the sterilization time. The AUTOMATIC MODELS are equipped with a digital microprocessor-based control (5),(11) which operates with a electronic type sensor (34) located near the steam exit, have too an automatic electrovalve (41) for the steam release and the capacity to program and store four different sterilization cycles with automatic control of temperature and time.

### ANALOGIC STERILIZER OPERATION

Be sure that control knob (11) is in OFF position before plug-in the equipment and before follow the next steps.

- a.- Close the drain valve (8)
- b.- Fill up the chamber with demineralized water up to the rack base level.
- c.- Put the rack (2) with the material inside the chamber.
- d.- Turn the lid until the arm (23) be inside the arm guide (24), check that lid gasket be inside his place all around, then turn the wheel (16) for tight well the lid.
- e.- Open completely the relief valve (20).
- f.- Site the control knob (11) in HI position ( maximum heating) for heat the water fast and begin

NOTE: Is important for a reliable sterilization that the air inside the chamber be removed from there. For this purpose the relief valve (20) must stay open at the beginning of each sterilization cycle and let the steam goes out for approximately 5 minutes before you close it. With this we are sure that the steam in the chamber is saturated, and temperature-pressure relation is the right.

**g.-** After 5 minutes of steam dragging close the relief valve, then the equipment will begin to pressurize. It is important to check the pressure in the gauge, when it's the desired value turn the knob toward a low heating position "LO", fix the knob in "2" position.

**h.-** Wait a moment, if the pressure is not the desired make the necessary adjustments in the knob. ("9" higher heating potency", "1" lower heating potency). Count yourself the sterilization time and then turn OFF the Sterilizer with the control knob.

**i.-** Before to open the Sterilizer is very important be sure that there is no pressure inside the chamber. Using a security glove open the relief valve to release the steam of the chamber; the pressure gauge will gradually descend to zero, and no more steam must come out of the relief valve. Let the sterilizer cool down a moment before opening.

IMPORTANT NOTE: This Sterilizer it is not an automatic, for this reason you must have permanent attention to the pressure gauge. In case that knob is fixed in a higher heating position than the

## AUTOMATIC STERILIZER OPERATION

Be sure that pilot switch (35) is in OFF position before follow the next steps.

Follow the same steps of the analog sterilizer (above cited) from **a)** to **d)**.

**e.-** Close the relief valve (20) and turn On the pilot switch (35).

The Sterilizer has an auto-check program that let you know from the most common failures of the system, therefore, the user can be sure that the Sterilizer is working in optimal conditions. Every time you turn on the Sterilizer, the next sequence is done: the display present for two seconds the number "8888" at the same time a beep is activated. After this, the displays turn off and the system starts the internal test and check too that the electrovalve is working. If all the features are right then the Sterilizer will present in the display in a intermittent way the last used program with temperature and time. If any of the features fail, the display will present the message "ErrX". Where the X could be a number from 1 to 6 depending on the error detected and the audible alarm will turn on. Press the increase button (  $\Delta$  ) to re-start, if the error persists, you should contact Felisa or any authorized distributor.

Failure codes:

Sensor failure	Err (1 and 5)
Heating Element failure	Err (6)

#### **f.- Sterilizing cycles Programing.**

The Sterilizer has a simple menu that allows you to configure up to 4 sterilizing programs (P0, P1, P2 and P3) each one with his temperature, time and evacuation type.

Press the menu button ( ☐ ), the display will present "ProG" in the upper display and the program number in the lower display, with the increase (  $\Delta$  ) or decrease (  $\nabla$  ) buttons you can select the desired program. Once you have selected the program only press the "start" button for begin that sterilizing cycle, or press again the menu button to modify his parameters as is described below:

##### **Setting the working temperature (SP)**

The display will show the last temperature for this program, with the increase or decrease buttons you can change to the desired working temperature, once you change it press the menu button

##### **Setting the overshoot alarm**

When you finish the configuration of the working temperature, the display will present the message "ALHi" and the last value of the alarm. This temperature will always be at least 3 ° C higher than the Set Point temperature. With the increase or decrease buttons you can modify the overshoot alarm value, once you have the desired value press the menu button to record and step ahead in the configuration.

##### **Sterilization time**

The upper display will present the message "timE" and the lower display will show the last time programmed, with the increase or the decrease buttons you can select the time wanted from 5 to 99 minutes. Once you made your selection press the menu button to record and go ahead in the configuration.

Note: It is important to know that the time start running once the Sterilizer has reach the programmed temperature (SP).

##### **Evacuation type**

The Sterilizer allow you to select the speed of the steam release in the chamber after the sterilization. Fast: The relief valve (20) is open permanently to evacuate the chamber in one minute approximately. Slow: The relief valve is open in short periods depressurizing the chamber in 10 minutes approximately. When you get in to this part of the menu, the upper display will present the message "EUAC" and the lower display the last selection ("FS" fast, or "SL" slow). With the increase or decrease buttons you can select the type of evacuation and then press the menu button.

##### **Audible alarm**

The Sterilizer has an audible alarm that the user can activate, when is activated will sound every time that a failure exist in the operation of the system. The display will present the message "Al" an then the status of the alarm ("ON" or "OFF"), with the increase or decrease buttons you can select your choice, once you made it press the menu button to record and finish the program setting.

#### **g.- Sterilization Process**

Once you have selected the sterilization program and it ´s tilting in the displays, to initiate the process push the "START" button, the Sterilizer will present in the upper display the actual

Is important remember that temperature sensor is located in the steam exit, not inside the water, for this reason you will see only a few rise in temperature until steam is generated and can heat the sensor.

The Sterilizer will be generating steam, this steam will drag the cold air out of chamber until the system detect no more increase in temperature, then will close de relief valve and a beep can be heard . The equipment keep on heating and pressurizing the chamber until it reach the programmed temperature (SP).

When the SP temperature is reached the "Sterilizing" LED will light and the lower display shows the countdown of programed time. After the sterilizing time is over a beep will sound, the equipment stops the heating, give a time for cool a few the temperature showing a line passing through the display. When the display presents "----" the equipment starts the chamber evacuation according the program.

**Is IMPORTANT do not open the Sterilizer lid until be sure that all the steam be out of the chamber and check the pressure gauge read zero.** When both displays flash showing the program again is the signal that the Sterilizer is ready for a new cycle.

#### **MAINTENANCE.**

Like any other manufactured product, some parts of the Sterilizer could be damaged after using for a long time, to replace them always use original factory parts, a list of parts is included in this instruction manual. All the parts can be ordered from any authorized Distributor or directly with Us.

#### **IMPORTANT NOTES.**

- To carry out any maintenance work disconnect before the Sterilizer from its energy source.
- Always check to have the proper water level in the chamber before start a new cycle, otherwise the heating elements will burn.
- Do not change the position of the temperature sensor.
- Connect the Sterilizer always to a ground- connected socket, otherwise it might cause electric shocks.
- Variations in voltage can damage the electronic components.
- If you need to stop the sterilization process for any reason and chamber is pressurized press the "START" button and the system will stop the process and will start to evacuate the steam as

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